

Surround Recommendations:

N-7 Paint Formula

Recommended for interior walls

Behr Paint available at Home Depot

Base – **Premium Plus Interior Flat (1500A)** pastel base

Formula ColorantOZ4896

B	Lamp black	0	15	1
I	Brown oxide	0	3	1
T	Medium Yellow	0	1	1

Recommended for the metal walls of a SpectraLight booth

Polyurethane

Sherwin Williams

Tel: 201-933-3800

Formula F63VXA3082-4350

Foam for viewing surface:

Uniroyal Engineering Products

2 North Tamiami Trail, Suite 900

Sarasota, FL 34236

Contact located in Stoughton, WI

800-873-8800 (Sue Deneen Ext 222)

X-Rite recommends servicing and re-certifying the Macbeth Lighting SpectraLight III once each year. This will ensure that your light sources meet manufacturers recommended specifications.

To properly repair and certify X-Rite instrumentation, we highly recommend the use of factory trained service engineers. Use original equipment manufactured parts for optimum performance.

Visit xrite.com for more information about X-Rite products. X-Rite customers worldwide may also call Customer Service at 800-248-9748.