

X-Rite Ci4200/Ci4200UV Service Programs



Optimize Your Spectrophotometer Investment

X-Rite's Ci4200 and Ci4200UV spectrophotometers are manufactured to rigorous quality standards. These standards are backed by a one-year standard warranty* that includes comprehensive global services and superior phone and web support. Now you can be assured that your device is always working at optimum performance when you choose one of the Ci4200/Ci4200UV Service programs.

Limit your operational down time, eliminate hassle, and ensure you are completely satisfied with your X-Rite products, year after year. Ideal for those requiring annual ISO or other certification and protection against unforeseen incidents.

**unless different regulations or custom contracts apply*



Choose Your Ci4200/Ci4200UV Protection Plan

Preventive Maintenance and Certification Contract

The Maintenance Contract provides for a one-time certification for the X-Rite Ci4200 and Ci4200UV. It is recommended to be recertified annually. Where applicable or necessary, X-Rite will provide printed documentation and/or a dated certification sticker.

Entitlements

One annual certification per contract. This certification meets the requirements of ISO 9001 and 17025. An X-Rite certified technician will annually:

- Verify and document the device's incoming performance
- Inspect and clean the device
- Remove cosmetic defects
- Exchange lamp if necessary
- Test measurement accuracy
- Verify repeatability & consistency
- Certify the device

Limitations

This plan does not cover unit repair. If a repair is necessary, it will be billed at either a flat rate or for labor and parts at the time of service. Consumables (other than lamp) and accessories are excluded. Damage caused to devices by abuse or misuse is not covered. **Ink, paint, or other contaminants in the sphere are not covered.**

Flat Fee Repair

The Flat Fee Repair plan provides one-time service for X-Rite Ci4200 and Ci4200UV devices out of warranty and not already under any X-Rite service contract. Where applicable or necessary X-Rite will certify your device and provide printed documentation and/or a dated certification sticker.

Entitlements

- Service work billed at flat rate fee
- Includes repair and replacement of parts identified as failed
- 90 day warranty on parts repaired/replaced

Full Service Contract

The Service Contract is available to X-Rite Ci4200 and Ci4200UV owners from the time of purchase until the expiration of the original product warranty. Contracts are available for one, two or three years.

Entitlements

- One annual certification per contract
- Labor and parts for repairs (excluding consumables and accessories)
- Repair warranties are valid for the length of the contract
- Outbound shipping included
- Loaner units available with advanced delivery and outbound shipping included
- Extends support via phone & email for the length of the contract
- Access to online support center
- Provides price protection for the length of your contract

Limitations

Consumables (other than lamp) and accessories are excluded. Damage caused to device by abuse or misuse is not covered. **Ink, paint, or other contaminants in the sphere are not covered.** Inbound shipping not included. Loaner units are offered on a first come, first served basis, and are subject to availability.

Units can be placed under contract while under warranty or if they have been certified or repaired within the last 30 days, if not under warranty a chargeable inspection/certification is required.

Bundles with Service Offerings

X-Rite offers bundled packages that include the Ci4200/Ci4200UV hardware along with NetProfiler and X-Rite Quality Control software. Please contact your local X-Rite sales representative for details.

Ci4200/Ci4200UV Service Options Summary

Value Proposition	Feature	Preventive Maintenance and Certification Contract	Full Service Contract
		PM4200-01	SC4200-01 SC4200-02 SC4200-03
Protect your investment	Extend original HW warranty for 1 year (parts and labor)	X	✓
	Price protection for the length of your contract	✓	✓
	Annual Certification of one unit	✓	✓
	Parts and labor included for repairs	X	✓
	Loaner Units with advanced delivery	X	✓
	One-way outbound shipping costs included	X	✓
Optimize and improve quality control	Includes phone and email technical Support	X	✓
	Access to online support center	✓	✓

X-Rite is the world leader in providing solutions for professionals involved with all aspects of color control. We lead the industry in offering service options to ensure uninterrupted performance for your Spectrophotometer. Training and educational resources are available globally and on line for both new and experienced users to optimize their color control capabilities.

X-Rite Training and Education

X-Rite offers a series of educational programs that consists of classroom, eLearning, regional seminars, personal one-to-one and on-site programs. These programs are designed to enhance your knowledge of color management technology and workflows.

For additional information on support programs or training and education, visit www.xrite.com/services or call your local service center.

Service Locations

Americas:

X-Rite, Inc.

Phone: 800 248 9748 or
+1 616 803 2100

Europe, Middle East & Africa:

X-Rite Europe GmbH

Phone: 00 800 700 300 01 or
+41 44 842 24 00

X-Rite - Munich, Germany

Phone: 00 800 700 300 01

X-Rite Italy S.r.l.

Phone: 00 800 700 300 01

X-Rite United Kingdom

Phone: 00 800 700 300 01

Asia Pacific:

X-Rite Asia Pacific Limited

Phone: +852 2568 6283

X-Rite China Service Center

Phone: +86 21 644 81155

X-Rite Japan

Phone: +813 6825 1641

X-RITE WORLD HEADQUARTERS

Grand Rapids, Michigan USA • (800) 248-9748 • +1 616 803-2100 • xrite.com
© 2013, X-Rite, Incorporated. All rights reserved. L10-480 (06/13)

